

Over-The-Wire Stoma Measuring Device


- To more accurately identify stoma tract length for better fit, comfort and safety
- Radiopaque allows fluoroscopic visualization
(Only contained in MIC-KEY G and MIC-KEY* J/TJ Introducer Kits)*

Serial Dilator with Peel-Away Sheath

- Smooth and tapered to minimize tissue trauma
- Incorporated Peel-away Sheath to enhance procedural efficiency
- One dilator in every kit to facilitate the correct dilation of the stoma

SAF-T-PEXY* T-fasteners

- Pre-loaded
 - No extra manipulation
- User friendly
- Resorbable 3/0 Biosyn™¹ suture:
 - No return to hospital needed for removal
 - Minimizes the risk of infection
- Locking polyurethane suture lock
 - Secures the stomach to the anterior abdominal wall



Hemostat

- To close the T-fastener suture lock

Safety Scalpel

- To create a small skin incision through the subcutaneous tissue and the fascia

Over-The-Wire Seeking Catheter

- Aids in manipulating the wire through the pylorus and advancing it past the Ligament of Treitz
(Only contained in MIC J/TJ and MIC-KEY* J/TJ Introducer Kits)*

Safety Introducer Needle

- To puncture the stomach to allow the j-tip guidewire to be advanced into the stomach
- Safety collar to prevent inadvertent needle sticks

Guidewire

- To maintain access and facilitate dilation
- J-tip to prevent perforation of the opposite stomach wall